

| 70Xi Series Lab Analyzers | | | | | Manual Pipette | | Automatic "Single Shot" Injection | | Autosampler System | |
|----------------------------|---------------------------|---------------------------|---|--------------------------------|------------------------------|---|-----------------------------------|----------------------------|---|-----------------------------|
| | | | | | Model | Available Options | Model | Available Options | Model | Available Options |
| Freeze Point ASTM D5972 | Cloud Point ASTM D5773 | Pour Point ASTM D5949 | | | PSA-70Xi | Antifreeze Freeze Point (ASTM D6660) Solid Point Crude Oil - Pour Point Black Oil - Pour Point | PSA-70XiSS (Fuels Only) | Solid Point | PSA-70XiAS (Fuels Only) ASL-70Xi (Fuels & Lubes) | Solid Point |
| Freeze Point ASTM D5972 | Cloud Point ASTM D5773 | | | | FCA-70Xi | Antifreeze Freeze Point (ASTM D6660) | FCA-70XiSS (Fuels Only) | | | FCA-70XiAS (Fuels Only) |
| Freeze Point ASTM D5972 | | | Pour Point ASTM D5949 | | PFA-70Xi | Antifreeze Freeze Point (ASTM D6660) Solid Point Crude Oil - Pour Point Black Oil - Pour Point | PFA-70XiSS (Fuels Only) | Solid Point | PFA-70XiAS (Fuels Only) | Solid Point |
| | | Cloud Point ASTM D5773 | Pour Point ASTM D5949 | | PCA-70Xi | Solid Point Crude Oil - Pour Point Black Oil - Pour Point | PCA-70XiSS (Fuels Only) | Solid Point | PCA-70XiAS (Fuels Only) APC-70Xi (Lubes Only) | Solid Point |
| Freeze Point ASTM D5972 | | | | | FPA-70Xi | Antifreeze Freeze Point (ASTM D6660) | FPA-70XiSS (Fuels Only) | | | FPA-70XiAS (Fuels Only) |
| | | Cloud Point ASTM D5773 | | | CPA-70Xi | | | CPA-70XiSS (Fuels Only) | CPA-70XiAS (Fuels Only) | |
| | | Pour Point ASTM D5949 | | | PPA-70Xi | Solid Point Crude Oil - Pour Point Black Oil - Pour Point | | | PPA-70XiAS (Fuels Only) APA-70Xi (Lubes Only) | Solid Point |
| Freeze Point ASTM D5972 | Cloud Point ASTM D5773 | Pour Point ASTM D5949 | Viscosity @ -20°C & -40°C ASTM D7945 | Density @ 15°C ASTM D4052 | | | JFA-70XiPC (Fuels Only) | Solid Point | JFA-70XiPC-AS (Fuels Only) | Solid Point |
| Freeze Point ASTM D5972 | Cloud Point ASTM D5773 | Pour Point ASTM D5949 | Viscosity @ 40°C ASTM D7945 | Density @ 15°C ASTM D4052 | | | DFA-70XiFP (Fuels Only) | Solid Point | DFA-70XiFP-AS (Fuels Only) | Solid Point |
| Freeze Point ASTM D5972 | | | Viscosity @ -20°C & -40°C ASTM D7945 | Density @ 15°C ASTM D4052 | | | JFA-70Xi (Fuels Only) | | | JFA-70Xi-AS (Fuels Only) |
| | | Cloud Point ASTM D5773 | Pour Point ASTM D5949 | Viscosity @ 40°C ASTM D7945 | Density @ 15°C ASTM D4052 | | | DFA-70Xi (Fuels Only) | Solid Point | DFA-70Xi-AS (Fuels Only) |
| Custom Application | | | | | GPA-70Xi | Freeze Point Cloud Point Pour Point Antifreeze Freeze Point Solid Point Gel Point Titer & Melting Point | | | | |

Options for all 70Xi Series: Phase Temperature Control Software (PTCS)
LLP-70A Label Printer
BSC-70X Barcode Scanner

| Portable Cloud Point Analyzer | | Manual Syringe Injection | |
|-------------------------------|--|--------------------------|--|
| | | Model | Available Options |
| Cloud Point ASTM D7397 | | | BAT-T30 Rechargeable Battery Pack DCA-T30 DC Power Adapter (with cigarette plug) |

| Wax Appearance Temperature Analyzer | | Manual Syringe Injection | |
|---|--|------------------------------|--|
| | | Model | Available Options |
| Wax Appearance Temperature (WAT) or Cloud Point and Wax Disappearance Temperature (WDT) | | WAT-70Xi (Crude Oil Only) | PTCS LLP-70A Label Printer BSC-70X Barcode Scanner |